BL21-Star(DE3)pLysS 感受态细胞

BL21-Star(DE3)pLysS Chemically Competent Cell

Cat.NO. ZC1208

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Sample Induction Protocol (for reference only)

1. Inoculate a single colony from a freshly streaked plate into 3ml of LB medium containing the appropriate antibiotic for the plasmid and host strain.
2. Incubate with shaking at 200rpm at 37 $^{\circ}$ C overnight.
3. Inoculate 50ml of LB medium containing the appropriate antibiotic with 0.5ml of the overnight culture prepared in step 2 (use the 500ml triangular flask as the container would be better).
4. Incubate with shaking at 150rpm at 37 $^{\circ}$ C until the OD₆₀₀ reaches 0.5-0.8.
5. (Optional) Pipet 1ml of the cultures into clean microcentrifuge tubes and place the tubes on ice until needed for gel analysis or storage at -20 $^{\circ}$ C. These will serve as the non-induced control samples.
6. Add IPTG to a final concentration of 1mM. Optimal time for induction of the target protein may vary from 2-16hours, depending on the protein.
7. Incubate with shaking at 120 rpm at 37 $^{\circ}$ C for 2-4hours. To determine the optimal time for induction of the target protein, it is recommended that a time course experiment be performed varying the induction from 2-16 hours.
8. Place the culture on ice for 10 minutes. Harvest cells by centrifugation at 5,000 \times g for 10minutes at 4 $^{\circ}$ C.
9. Remove the supernatant and store the cell pellet at -20 $^{\circ}$ C (storage at lower temperatures is also acceptable).

IPTG 配制:

Prepare a 1M solution of IPTG (Isopropyl- β -D-thiogalactoside; Isopropyl- β -D-thiogalactopyranoside) by dissolving 2.38g of IPTG in dd water and adjust the final volume to 10ml. Filter sterilize before use.

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Chloramphenicol 34mg/ml in ethanol. Store at -20 $^{\circ}$ C. Use at 34 μ g/ml.