# ArcticExpress(DE3) pRARE2 感受态细胞

## ArcticExpress(DE3) pRARE2 Chemically Competent Cell

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### Sample Induction Protocol (for reference only)

1. Inoculate a single colony from a freshly streaked plate into 3ml of LB medium containing the appropriate antibiotic for the plasmid and host strain.

2. Incubate with shaking at 200rpm at  $37^{\circ}\text{C}$  overnight.

3. Inoculate 50ml of LB medium containing the appropriate antibiotic with 0.5ml of the overnight culture prepared in step 2 (use the 500ml triangular flask as the container would be better).

4. Incubate with shaking at 150rpm at  $37^{\circ}\text{C}$  until the  $\text{OD}_{600}$  reaches 0.5-0.8. (0.6 recommended; about 2.5h).

5. (Optional) Pipet 1ml of the cultures into clean microcentrifuge tubes and place the tubes on ice until needed for gel analysis or storage at  $-20^{\circ}\text{C}$ . These will serve as the non-induced control samples.

6. Add IPTG to a final concentration of 1mM. Optimal time for induction of the target protein may vary from 2-16hours, depending on the protein.

7. Incubate with shaking at 120rpm at  $37^{\circ}\text{C}$  for 2-4hours. To determine the optimal time for induction of the target protein, it is recommended that a time course experiment be performed varying the induction from 2-16hours.

8. Place the culture on ice for 10minutes. Harvest cells by centrifugation at  $5,000\times g$  for 10minutes at  $4^{\circ}\text{C}$ .

9. Remove the supernatant and store the cell pellet at  $-20^{\circ}\text{C}$  (storage at lower temperatures is also acceptable).

### IPTG 配制：

Prepare a 1M solution of IPTG (Isopropyl- $\beta$ -D-thiogalactoside; Isopropyl- $\beta$ -D-thiogalactopyranoside) by dissolving 2.38g of IPTG in dd water and adjust the final volume to 10ml. Filter sterilize before use.